

AMRAD™



AMRAD ADVANTAGE MEDICAL X-RAY SYSTEMS

Take a closer look

THE AMRAD ADVANTAGE ORTHOPAEDIC X-RAY

AMRAD'S™ fully integrated system is specifically designed to meet the demands of today's orthopaedic x-ray procedures. Amrad™ is known for innovative development and quality workmanship in radiographic systems and components. Amrad™ products help physicians provide top-quality health care for their patients every day.

ENGINEERING AND MANUFACTURING EXCELLENCE

Each component of the Amrad orthopaedic system — generator, table, tubestand and wallstand — has been designed for long term durability under heavy patient volume conditions. And, as always, our 5-year parts warranty attests to our skills and commitment.

The designs of the Amrad engineers have resulted in the manufacture and sale of over 20,000 high quality imaging systems since 1977. Quality standards are unsurpassed, yet the products remain the affordable choice in today's marketplace.

AMRAD'S WORLD CLASS SYSTEM COMPONENTS

- Deluxe four way float-top table with reciprocating bucky.
- Floor mounted transverse/rotational tubestand with tube trunnion capability.
- Conventional Wallstand or optional Deluxe Wallstand with *trayless*, self-centering feature.
- Choice of quality, high performance x-ray tube and collimator.
- Selection of dependable, user-friendly high frequency generators.
- All backed by a 5-year parts warranty.

DIGITAL IMAGING READY

AmradAdvantage radiographic equipment is completely compatible with most digital imaging systems. Digital imaging technology may be added at any time. Our focus on the future secures your investment today and guarantees your investment value tomorrow.



SYSTEM



FMT system with optional elevating table and Deluxe Wallstand.

AMRAD

THE AMRAD™ FLOOR MOUNTED TUBE STAND

Amrad's™ floor-rail rotational tubestand offers easy positioning for the most difficult views. The extended vertical travel of the tube to the floor allows for weight bearing ankle and knee views. Column rotation and tube trunnion capabilities facilitate cross-table exams, while 8" of transverse tube travel allow effortless extremity procedures. This smooth handling tubestand is extremely operator friendly. An operator control handle provides easy operation of all tubestand locks. Our tubestand is designed to accommodate almost any size facility and provides a stable, heavy-duty support structure. The simple floor-rail mounting reduces installation costs.



THE AMRAD DELUXE WALLSTAND

This unique wallstand features *trayless* cassette centering to reduce procedure time. Merely insert the cassette, depress the handle and the cassette automatically centers. Vertical travel of the wallstand reaches the floor for weight bearing limb procedures and travels upward for chest procedures on the tallest of patients. Convenient electric locks secure the cassette cabinet at any elevation. This time-saving wallstand is offered in right or left hand load configurations with optional automatic exposure control chambers.



THE AMRAD RADIOGRAPHIC TABLE

A versatile, high performance radiographic table designed for all general radiographic applications. Exceptionally smooth and quiet vertical travel of the optional elevating table feature provides effortless transfer of stretcher, wheelchair, casted or geriatric patients and accommodates extreme patient sizes. The four-way floating table top is perfectly flat and designed for maximum ease during patient positioning as well as performance of extremity procedures.

The unique design of the table pedestal minimizes joints and cavities to facilitate disinfection. Our effective foot pedal design for table top locks and elevator motion enables quick activation by the operator while minimizing the possibility of patient activation. A high quality reciprocating bucky and grid complete the table. Amrad's™ elevating table is well-suited to the daily demands of orthopaedic patient procedures.



WORLDWIDE SUPPORT

Amrad™ radiographic products are backed by an exclusive, worldwide dealer network for support and service. This exceptional group of professionals is dedicated to providing you with unequalled products and services.